

# Development of a Portable PEMF Therapy Device



**Department of Electronics & Communication Engineering  
National Institute of Technology Patna  
Bihar-800005, India**

Presented by:

Tushar pradhan

Anurag sharma

Tharak Amruth Raj



# Contents

- Introduction
- Working Principle of PEMF Therapy
- Proposed Device Overview
- Block Diagram
- Circuit Diagram
- Component List for PEMF Device
- Working Process of the PEMF Device
- Expected Results & Benefits
- Conclusion & Future Scope



# Introduction

Arthritis causes chronic pain, limiting mobility and quality of life. Pulsed Electromagnetic Field (PEMF) therapy is a non-invasive method that helps reduce pain and inflammation.

## **Purpose:**

- Develop a **portable, compact, and user-friendly** PEMF device for arthritis pain relief.
- Provide an **affordable and accessible** solution for daily pain management.

## **Objective:**

- Design a **lightweight and adjustable** PEMF device.
- Ensure **ease of use** for patients of all ages.
- Offer **effective pain relief** without medication or side effects.

# Working Principle of PEMF Therapy

## 1. Electromagnetic Field Generation

- The device generates pulsed electromagnetic fields (PEMF) that penetrate deep into body tissues.
- These fields induce small electrical currents within cells, stimulating biochemical and physiological processes.

## 2. Cellular Stimulation & Healing


- Enhances ion exchange: Restores cellular balance, improving energy production (ATP synthesis).
- Increases blood flow: Promotes oxygen delivery and nutrient absorption for faster healing.
- Reduces inflammation: Modulates inflammatory responses, easing joint stiffness and swelling.
- Accelerates tissue regeneration: Encourages the repair of damaged cartilage, muscles, and bones.

## 3. Scientific Validation

- NASA research confirms 1-4 Gauss as the optimal intensity for cell regeneration and bone healing.
- Studies show PEMF therapy reduces pain by blocking pain signals and increasing endorphin production.

## 4. Non-Invasive & Drug-Free Therapy

- Offers safe and painless arthritis relief without side effects.
- Can be used daily for long-term pain management and mobility improvement.



# Proposed Device Overview

## **1. Controlled PEMF Signal Generation**

- Produces precise and adjustable pulsed electromagnetic fields for effective pain relief.
- Ensures optimal intensity (1-4 Gauss) for safe and efficient therapy.

## **2. Portable & User-Friendly Design**


- Compact, lightweight, and easy to carry for home or travel use.
- Simple interface with one-touch operation for convenience.

## **3. Adjustable Intensity & Frequency**

- Customizable settings to suit individual pain levels and treatment needs.
- Offers multiple frequency options for targeted therapy.

## **4. Medical Applications**


- Helps manage arthritis, joint pain, muscle stiffness, and inflammation.
- Potential use in post-surgical recovery, bone healing, and sports injuries.



# Block Diagram:-


Power source

- Signal generator (PC, mobile app, or recorded audio)
- Signal amplifier
- Output coil for electromagnetic radiation
- Optional feedback mechanism (magnetometer or ammeter)



## Circuit Diagram :-

Visual representation of the circuit with components and their connections.





# Component List for PEMF Device :-

## **Active Components :-**

- IRF3205 MOSFET x3 → Handles high current switching.
- 2N2222 NPN Transistor x1 → Controls MOSFET switching.
- Rotary Potentiometer (10K $\Omega$ ) x1 → Adjusts intensity.
- 1N4007 Diode x1 → Protects against voltage spikes.

## **Passive Components :-**


### Resistors:

- 100 $\Omega$  x3 → Limits gate current for MOSFETs.
- 2.21k $\Omega$  x1 → Biasing resistor for transistor.
- 1k $\Omega$  x1 → Controls base current for transistor.

Coil (Radiation Coil) → Generates pulsed electromagnetic field.

### Other Essentials:

- Power Supply → 12-24V DC (Battery/Adapter).
- Soldering Iron & Solder Wire → For circuit assembly.
- PCB or Breadboard → For circuit prototyping.
- Connecting Wires & Casing (PVC Board) → For connections and enclosure.





# Working Process of the PEMF Device :-

## **1.Signal Input**

- The input signal is generated from a PC, mobile device, or an external PEMF generator.
- This signal determines the pulse frequency and intensity.

## **2.Signal Amplification**

- The 2N2222 transistor acts as a switch, controlling the MOSFETs.
- IRF3205 MOSFETs amplify the signal to deliver sufficient current to the coil.
- The potentiometer allows adjustment of intensity as per user needs.

## **3.Coil Emission**

- The amplified signal is sent to the radiation coil.
- The coil converts this electrical energy into a pulsed electromagnetic field (PEMF).

## **4.Therapy Application**

- The PEMF waves penetrate body tissues, stimulating cellular repair and reducing pain.
- This non-invasive therapy helps reduce inflammation, improve blood circulation, and enhance healing.

# Expected Results & Benefits :-



## ✓ Portable & Affordable

- Compact design makes it easy to use anywhere.
- Cost-effective alternative to expensive commercial PEMF devices.

## ✓ Customizable Therapy Settings

- Adjustable intensity and frequency to suit individual needs.
- User-friendly controls for personalized pain relief.

## ✓ Effective Pain Management

- Helps reduce arthritis pain, inflammation, and stiffness.
- Promotes faster healing and improved mobility without medication.

# Conclusion & Future Scope :-

## ✓ Conclusion:

- The proposed portable PEMF device offers an affordable, non-invasive solution for arthritis pain relief.
- It provides customizable therapy settings, ensuring effective and user-friendly treatment.
- The device has the potential to improve quality of life for arthritis patients by reducing pain and inflammation.

## ✓ Future Scope:

- Adding a Feedback Meter: Real-time monitoring of therapy intensity and effectiveness.
- Smartphone App Integration: Wireless control via a mobile app for convenience.
- Enhanced Coil Designs: Improved efficiency and better penetration for optimized therapy results.

**THANK YOU**

